Nanoscribe Photonic Professional GT 3D Lithography System



The Nanoscribe Photonic Professional GT system utilizes two-photon polymerization for 3D micro printing and maskless lithography. Printing of 3D structures can be achieved using conventional direct laser writing (DLW) or Nanoscribe's dip-in laser lithography (DiLL). The system is equipped with a piezo for printing in high resolution arbitrary 3D trajectories and a galvanic mirror for fast layer-by-layer printing. Features 200 nm in size to structures of several millimeters can be printed. Significant reduction of printing time can be achieved with the "scaffold and shell printing" method with post development UV curing. The system accepts STL files which can be created using CAD software. The files can be converted to Nanoscribe's GWL code using its Describe software. Users may also opt to create structures by writing GWL code directly. Nanowrite software allows for easy operation of tool with AxioVision software allowing for live monitoring of the process. The system can be used for a broad range of applications on the nano-, micro-, and mesoscale including; from optics, plasmonics, microfluidics, sensors, micro-robotics, biomimetics, to life sciences, and more.

Technical Specifications:

- Objectives 63X immersion, 25X immersion, 10X Immersion and 20X (for air)
- Laser Class 3 120 mW, 780 nm, 150 fs, 80MHz fiber laser
- 100 x 100 mm2 motorized XY stage range
- 300 x 300 x 300 μm piezo range for 25X and 1mm diameter for 10X
- Sample holders: DiLL (dip-in laser lithography) 1 X 1-inch glass slides, users can bring their own substrates
- Nanoscribe's Proprietary Photoresists: IP Dip, IP-L 780, IP-S, IP-Q, IP-Visio
- Dymax BlueWave 200 curing system available for curing of structures post-writing



Parameter Edits

You will need to edit the parameters to be used for your print job. The editing is done in Nanoscribe's Describe software and is converted to GWL code upon saving the file. A detailed description of each parameter is given in the appendix. This manual only covers Dip in Lithography (DILL). Other forms of printing are cover in supplementary manuals that the user can request.

DILL printing uses the 25x objective along with IP-S. It can also be done with the 63x objective along with IP-DIP. The objective moves away from the substrate during printing so ITO covered glass or other substrates can be used.

- 1. Open your .STL file in the Describe software Descri
- You will see the following page come up. From the Input dropdown menu, select the 25x objective. Once done, click Slice at the top of the page. If you wish, change the size of the structure using Scaling. The Lock Aspect Ratio option can be used to scale the structure in all 3 dimensions simultaneously.

Mesh statistics Number of triangles	Model Slice >							
40334 Dimensions X : 100 μm Y : 77.697 μm Z : 82.453 μm Volume 7.12×10^-5 mm ³ Display options Disable preview	Input beefyshark-medium.stl Browse Last session parameters Default parameters IP-Dip 63x Solid (Fused Silica) IP-S 25x Solid (ITO) IP-S 25x Shell (ITO)							
Options	Orientation							
	Scaling							
a	Axis Absolute (µm) Relative							
N	X 100 \$ 1.40446 \$							
2	Y - 77.697 🗘 1.40446 🗘							
	Z 82.453 \$ 1.40446 \$							
10 γ (μ), μα 10 10 10 10 10 10 10 10 10 10	Lock aspect ratio							

3. In the **Slice** window, select the **Slicing Mode** you want. Also, enter the maximum and minimum **Slicing Distances** you want. Click **Slice** at the bottom to save any changes. Once done, click **Fill** at the top.





Slicing stats	< Model	Slice	Fill >
Number of slices 405 Number of points 53047	Slicing parameters	ters	
Display options	◯ Fixed (●)	Adaptive	
Disable preview	Slicing distar	nce	0.5 🖨
Options	Minimum slici	ing distance	0.2 🖨
2000	Slope evalua	tion	
10 10 10 10 10 10 10 10 10 10	Average		Max
10 10	Surface norm	nals	
	Smooth		Flat
2 (m) (n	Voxel aspect	ratio	6 ‡
	Simplificatio	n tolerance	0.05 🗘
	▼ Fix self-in	tersections	
Ψ (μt) 100 100 100 100 100 100 100 100 100 10		Slice	

4. In the **Fill** window, select the **Fill Mode** you want. Also, changes any other parameters you want. Click **Generate Shell** at the bottom to save any changes. Once done, click **Scaffold** at the top.

Display options Clipping	✓ Slice Fill	Scaffold $>$		
Axis © X Y Z	Fill mode Solid Shell & Scaffold			
 ✓ Show shell ✓ Show mesh overlay 	Shell parameters Hatching distance	0.5 🖨		
Options	Shell contour count	12 🖨		
101 m 101 m 101	Base slice count Hatching angle	2 🖨		
	0° 0 °	0 \$		
· · · · · · · · · · · · · · · · · · ·	Hatching angle offset			
	▽	• •		
20 V (a0 11 10 10 10	Generate shell			

UNIVERSITY OF ILLINOIS AT CHICAGO CORE FACILITY COLLEGE OF ENGINEERING In the Scaffold window, select the Scaffold Paramters you want. Changes any other parameters you want. Click Generate Scaffold at the bottom to save any changes. Once done, click Output at the top.



 In the Output window, select the Overwrite Existing. Changes any other parameters you want. Once done, click Save at the top. Close the Describe software.

2000 180	< Scaffold	Output Save >	Array
188 X (µn) 160 140	Output options - Base file name	beefyshark-mediu Browse	-Splitting -Splitting mode
120	Job file	beefyshark-medium_job.g	None Rectangular Hexagonal
80	Data file	beefyshark-medium_data.	X 70 ¢ Y 70 ¢ Z 70 ¢
6e			Block offset
60 40 Z (un)	Scan mode	O Piezo	X 142.5 ¢ Y 142.5 ¢ Z 0 ¢
Z (m) m		Calvo	Block width
28	Z-Axis	Piezo	X 81.6 μm Y 91 μm Z 71 μm
20		Microscope z-drive	Block shear angle 15 🗘
	Exposure	Different for shell & scaffold	Block overlap 2 🗘
200		© Constant	Layer overlap 1 🗘
180	Z-Direction	Opwards (+Z)	Block order
		O Downwards (-Z)	Avoid flying blocks
Y (µm) 1að	Hatch lines	● ↑↓ ① ↑↑	Group neighboring blocks
120			Backlash correction

UNIVERSITY OF ILLINOIS AT CHICAGO COLLEGE OF ENGINEERING In new window that comes up, you will see the GWL code that has been generated. Click the 3D Preview, then click Generate 3D Preview. A preview will be generated. You can click play at the bottom to watch an animation of the printing process.

💋 Ei	ffelTow	er (repariert)_	job.gwl -	DeScribe 2	.5.4	_		
File	Edit	3D Preview	Debug	Window	Help			
	EiffelTo	Genera	te 3D Pre	view	F	:5		Ŧ
1	% Fi	Regene	erate last 3	3D Preview	S	Shift+F5	· · · ·	-
3	% Sy Inve	Fit Visil	ble		F	7		
6	% Wr	Switch	to Standa	rd View		•		
7	Gal\ Cont	Render	ring Mode	•		•		
10 11	Galv	Save So	creenshot		F	-12		
12 13	% Sc	Quick-	Save Scree	enshot	C	Ctrl+F12		
14 15	XOFI YOFI	3D Pre	view Optio	ons				



UNIVERSITY OF ILLINOIS ATCHCAGG OF ENGINEERING

Sample Preparation and Writing

The sample holders are in the cabinet next to system. This section of the manual will cover how to use 25x objective for DILL with IP-S photoresist on ITO coated glass. The preparation for other resists for other printing modes is similar.

1. Open the cabinet and take out the sample holders, tape, and sample holder stand. Depending on the size of your sample, choose the holder appropriate for you.



Place the sample holder on the sample holder stand with the part labeled BOTTOM {RESIST} facing up towards you.



3. Place the sample in the appropriate shape. For ITO coated glass piece, make sure the ITO coated part of the sample is facing up toward you. Tape down the sample.





NOTE: To find the side with the ITO on the glass sample, use a multimeter to measure the resistance of the sample.



Pick up the photoresist bottle and wipe off the excess on the bottles edge. Place a drop of resist on the glass.
 You don't need a lot of resist. Let the resist spread out on the glass for about 5 minutes.



UNIVERSITY OF ILLINOIS ATCHCAGO CORE FACILITY COLLEGE OF ENGINEERING



5. Open Nanowrite software Nou will be asked if you want to calibrate the stage. Click Calibrate and allow

initialization to complete. Once done, the message **Initialization Done** will appear in the Nanowrite dialog box.

	Stage calibration	×
	Please remove any obstacles out of the movement a of the stage and make sure that the microscope dis indicates "lower z-limit reached" before clicking [Ca CAUTION: The movement of the stage cannot be st and might lead to serious injuries! Calibrate Disable stage	area play librate]. opped
NanoWrite 1.8.1	1	
System Idle	Camera Advanced Settings GWL mini script Submit	Messages ? Intralizing System Galvoscanner initialized in analog mode! ^ Initializing interlocks Initializing inlumination
Approach Sample		S08 VIS - 509 L LSM Duo T 80 / R 20 IR Objective Pos.: 5, i LCI Plan-Neofluar 25x/0.8 Imm Korr DIC M27 Initializing place
Load Job		Moving stage to default position E Cleaning up compilation cache Initialization done. Move to sample center Moving stage to default position
-Focal plane -	PerfectShape Interface Tilt Correction Skip Settings	Neady for exchange Noving stage to defailt position
Interface at (µm) 0.000 Find Interface	On/Off Transmission	X 000.004 0000.000 shape Y 000.001 0000.000 Invert Invert Z 000.006 0000.000 Step (µm) 1
Start Job	Reflection On/Off	
Abort	Progress 00:00:00 Estimated Time 00:00:00 Estimate Time	DILL

6. Click on the Exchange Holder tab in Nanowrite. Check the LCD monitor (Nanoscribe Control Display) on the table and make sure the message LOWER Z LIMIT REACHED is displayed. If not, turn the large knob on the side of the monitor anti-clockwise to display the message. Click OK on the Confirm Excange Holder dialog box.



UNIVERSITY OF ILLINOIS AT CHICAGO COLLEGE OF ENGINEERING

- Objective: LCI PInN 25x/0.8 W DICII Objec tives Re-Light Home Resolution: 0.42 µm Home flector path Optovar: 1x Tot.mag.: 250x HAL: 0.3V TL Reflector: Pos. -3-63x Oil position 20x sition Microscope **TL Illumination** Control TL Illumination Pos. 6 LSM Duo T 80 / R On w VIS: 50% Sideport: 50% left Z-Position: 0.00000 mm Settings Off On Automatic Pos. 5: LCI PInN 25x/0.8 W DICII Contrast Manager Make It Visible! Display Display BF LSM IR TROP
- 7. On the LCD monitor on the table, press Microscope and choose the objective lens you want to use for writing.

NOTE: SOMETIMES THE TURRENT MIGHT NOT ROTATE TO POSITION. IN THAT CASE, FEEL FREE TO ROTATE IT INTO POSITION BY HAND. BE GENTLE.

NOTE: THE 20X OBJECTIVE IS FOR AIR ONLY. DO NOT DIP IN PHOTORESIST OR OIL. PLEAESE TELL NCF STAFF IMMEDIATELY IF THIS HAPPENS SO THAT THE OBJECTIVE CAN BE CLEANED.

The Desktop should show this screen. PLEASE NOTE: DO NOT EXCHANGE, UNLOCK OR TOUCH SAMPLE HOLDER
 IF THIS SCREEN IS NOT DISPLAYED.





9. Hit "Open" on the interlock panel on the machine. Then lift the handle to open the loading area.



10. Gently open the locks on the piezo and lift up until it locks. Once locked, pull the pin on the right to make it completely vertical. MAKE SURE IT IS LOCKED. IF IT FALLS, IT CAN COST A LOT TO FIX.



11. Choose the lens you are using on the Nanoscribe control display and then screw the lens in the spot.



12. Once done, engage the locks of piezo and gently rest piezo on the locks for putting in the sample holder. **DO NOT FORCE THE HOLDER IN.**



- 13. Once sample is in, lift piezo and disengage locks of piezo to put piezo back and then lock it in place.
- 14. Close the lid and press gently until locked.
- 15. On the computer, choose the position your sample is in on the picture and click **OK**.



16. Click on Transmission/Illumination in Nanowrite.



17. Open the **AxioVision** software 4.8 and click on Live View.

	Axio	Vision	n LE															x
Fil	e I	Edit	View	Acq	uisition	Functions	Ann	otation	Measure	Eval	uate	Applicati	ons	Tools	Window	Help		
		0	R		y.	<u>,</u>	П,		¥		(1	۵	Η	Az			
d C	pen)	Imag	e Brow	/ser	Save As	Preview	ROI	Copy RO	I Cut	Сору	Paste	Live	Snap	Scale b	ar Text	Navigator	Workarea	Galler
Rea	dy																	

18. Click Approach Sample in Nanowrite. For ITO coated glass and IP-S resist, the software should find the interface automatically. Once found, you should see the message Interface Found, Strength: XX @ XX. You should also see the ripples in the Interface Finder window.



19. Click Load Job in Nanowrite and choose the desired job file. In your folder, click on the ABC_job.gwl file. You should see the message ABC_job.gwl Loaded in the message box.



20. If you have not found the interface of the sample and resist yet, click on **Find Interface** in **Nanowrite** find the interface.

NOTE: The microscope is equipped with an autofocus system to automatically focus the laser spot at the substrate surface. The autofocus system approaches a grating on the front focal plane of the microscope objective. The position of the image of the grating on the autofocus camera is an indicator of the position of the focal spot of the laser on the substrate. A detailed description of the autofocus system may be found in the supplied manual of the Definite Focus.

- 21. Click Start Job in Nanowrite to start printing. At any time, the job can be aborted by clicking Abort.
- 22. Once done. You will see the message Job Done. Click Exchange holder and ensure the message Lower Z Limit Reached appears.



23. Pull out sample holder.



DEVELOPING THE SAMPLE

1. Place a cleanroom wipe in the fume hood and label the wipe with your name, the name of the chemical (s), and the date. There is dedicated Nanoscribe glassware in the cabinet where the sample holders are located.



2. Place the sample in the special holder for glass slide and put it in SU-8 Developer. Once all the resist is gone, remove the sample and place it in IPA. The let the sample air dry. Do not use N2 gas to dry the sample. It could damage the structure(s) you just printed.



3. Put the chemicals in the waste bottle once done. Rinse the sample holder and glassware. Put glassware back into their box.

NOTE: You may bring your own sample holder for development. You may also bring your own glassware.

Cleaning Up and Shutting Down

- 1. Clean the sample holder with IPA.
- 2. Put the sampler holder and sample holder stage back in the cabinet. Put the resist back in the fridge.
- 3. Clean up the glass ware and rinse them with DI water.
- 4. Exit all software. Wait for the program to close completely (including all overlapping icons in the task bar)
- 5. Make sure the objective z-drive is completely lowered. If not, press on Load position on the touch panel and wait for the movement and the press Set work position.
- 6. Log out from the PC. **DO NOT SHUT DOWN.**